# MAYFIELD® 2 Skull Clamp

## **REF** (A3059)

### Instruction Manual

|                      |     |
|----------------------|-----|
| EN – English .....   | 3   |
| FR – Français .....  | 27  |
| IT – Italiano .....  | 51  |
| DE – Deutsch .....   | 75  |
| ES – Español.....    | 99  |
| NL – Nederlands..... | 123 |

This page intentionally left blank.

## MAYFIELD<sup>®</sup> 2 Skull Clamp

**REF** (A3059)

## Instruction Manual

### EN – ENGLISH

#### **United States, Canada, Asia, Pacific, Latin America**

USA 800-997-4868 ▪ 888-980-7742 fax

Outside USA 609-936-5400 ▪ 609-750-4259 fax

[integralife.com](http://integralife.com)

#### **Europe, Middle-East, Africa**

International: +33 (0)4 37 47 59 50 ▪ +33 (0)4 37 47 59 25 fax

Benelux: +32 (0)2 257 4130 ▪ +32 (0)2 253 2466 fax

France: +33 (0)4 37 47 59 10 ▪ +33 (0)4 37 47 59 29 fax

Germany: +49 (0) 2102 5535 6200 ▪ +49 (0)2 102 5536 636 fax

Switzerland: +41 (0)22 721 23 00 ▪ +41 (0)22 721 23 99 fax

United Kingdom: +44 (0)1 264 345 780 ▪ +44 (0)1 264 363 782 fax

[integralife.eu](http://integralife.eu)

#### **Manufacturer:**

Integra LifeSciences Corporation  
4900 Charlemer Drive, Building A  
Cincinnati, OH 45227, USA.

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69800 Saint Priest – France

Tel: 33 (0) 4 37 47 59 10

Meaning of Symbols used in this manual

 **CAUTION** 

Hazards which could result in equipment or property damage

 **WARNING!** 

Hazards which could result in severe personal injury or death.

Caution, consult accompanying documents

Product complies with the requirements of MDR 2017/745

Manufacturer

Authorized Representative in the European Community

**R<sub>x</sub> ONLY**

Caution: Federal (US) Law restricts this Device to sale by or on the order of a licensed healthcare practitioner

Catalog number

This device is not indicated for use in MR environment

Consult Instructions for Use

Medical Device

## Inspection

Always inspect instruments before and after use. If a component appears damaged and/or does not seem to function properly, do not use the device and immediately send the instrument to an authorized Integra repair center for evaluation, repair or replacement. Allow your Integra Representative to inspect this device a minimum of two times per year to assist you in its proper function.


Figure 1 A3059 Skull Clamp

1. 80 lb Torque Screw (454A1002)
2. Ratchet Extension
3. Rocker Arm
4. Release Button (both sides)
5. Index Locking Knob
6. Serial Number Label
7. Skull Clamp Base

The MAYFIELD® 2 Skull Clamp is designed for use with the following equipment:

- A3100 MAYFIELD 2 Base Unit
- A3101 MAYFIELD 2 Base Unit
- A3102 MAYFIELD 2 Base Unit
- A3018 MAYFIELD 2 Standard Swivel Adapter
- A3008 MAYFIELD 2 Tri-Star Swivel Adaptor
- ALL Styles of MAYFIELD Skull Pins (A1072, A1083, A1084, A1119, A1120, A1121, A1122, A1047, A1048)

Note: Use of MAYFIELD products and accessories in conjunction with other manufacturer’s stabilization equipment is not recommended.











|           |          |                                                                                     |                                                                                     |                                                                                       |                                                                                       |
|-----------|----------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Correct   | Figure A |    |    |    |    |
| Correct   | Figure B |   |   |   |   |
| Incorrect | Figure C |  |  |  |  |

Table 1 Skull Clamp Position Diagram

Note: Patient positioning and application of the skull clamp can be made with numerous variations. The patient’s final position should be acceptable to the surgeon, taking into account the patient’s anatomy, pressure points, neck flexibility and access to the surgical target.

In the chart above, the ovals representing the patient’s head have centerlines drawn on them. The diagrams show examples of proper and improper applications of the skull clamp. The direction of the single-pin should point directly through the center point created by the centerlines of the patient’s head with the pins on the two-pin side of the skull clamp being equidistant to the single-pin line of direction to obtain the highest degree of stability.

## **Indication For Use/Intended Purpose**

The MAYFIELD 2 Skull Clamp is placed on the patient's skull to hold their head and neck securely in a particular position when rigid fixation is desired. The Skull Clamp is indicated for use in open and percutaneous craniotomies as well as spinal surgery when rigid fixation is necessary.

The MAYFIELD 2 Skull Clamp is **NOT** Radiolucent.

### **WARNING!**

Failure to read and follow instructions furnished in this product insert may result in skull pin slippage and serious patient injury, such as lacerations, skull fracture, or even death.

The Skull Clamp must be applied along the center line of the patient's head with pins entering the skull perpendicularly (See Table 1). Failure to properly position the Skull Clamp on the patient's head could result in patient injury such as scalp laceration due to skull pin slippage. Failure to properly position patient and to fully secure all adjustment positions of the Skull Clamp or any support device may result in patient injury.

It is cautioned that only MAYFIELD products have been validated for use with other MAYFIELD products. Integra cannot advise whether other suppliers' products would function properly with the MAYFIELD line of products. Integra does not offer warranty protection when a non-MAYFIELD brand product is used with a MAYFIELD branded product.

Do not alter the design of the device in part or whole as serious patient injury could result.

Do not use Skull Clamps to exert excessive traction during a procedure. Patient injury could result.

Skull Clamps are not intended for use in or near the vicinity of a strong magnetic field (MRI).

### **WARNING!**

This device is NOT Radiolucent.

## **Intended Population**

MAYFIELD Skull Clamp fixation devices are not recommended for use on children under five (5) years of age. Extreme caution should be exercised in pediatric cases because of the thin skull.

## **Description**

The MAYFIELD 2 Skull Clamp is designed to provide rigid skeletal fixation or stabilization in conjunction with the Base Unit (A3100, A3101, A3102) or any other standard aluminum MAYFIELD Base Units (A2009, A2101, A2102). The Skull Clamp is designed to allow the surgeon freedom in positioning the fixation pins. Avoidance of critical areas of the skull is facilitated by a swiveling rocker arm that also rotates 360 degrees. To simplify patient re-positioning after skull pin impingements, the rocker may be rotated without adjustment of the impingement force imposed by the Torque Screw (Single-pin side of skull clamp).

The MAYFIELD 2 Skull Clamp (A3059) includes:

- Standard Adult Rocker Arm
- Standard Torque Screw (454A1002)
- Standard Ratchet Extension

## **Standard Torque Screw**

The torque screw can develop up to 80 lbs. (356 N) of force; graduated in 20 lb. (89 N) increments. The knob has eight facets. One facet is labeled with the “Caution” symbol.


Figure 2 Standard Torque Screw (454A1002)

## **Mounting of MAYFIELD Skull Pins**

Note: Use only MAYFIELD brand Skull Pins with this Skull Clamp.


Figure 3 Mounting Skull Pins

### **⚠ CAUTION! ⚠**

Carefully inspect reusable skull pins before each use. To maintain proper fixation, pins should be of equal length. Damaged (dull, bent, etc.) pins should be replaced. It is recommended that spare pins always be kept available.

### **⚠ WARNING! ⚠**

Never attempt to sharpen skull pins. Dull or bent pins should be replaced.

Note: When using MAYFIELD Skull Pins, follow Instructions For Use supplied with the product.

### **⚠ WARNING! ⚠**

Skull pins are not recommended for use on children under five (5) years of age and extreme caution should be exercised in pediatric cases because of the thin skull. MAYFIELD Child Skull Pins are recommended for young children over the age of five (5) years. It is at the discretion of the user, when the adult or child skull pins should be used.

### **⚠ CAUTION! ⚠**

Only MAYFIELD branded skull pins may be used with MAYFIELD Skull Clamps. Integra does not offer warranty protection when a Non-MAYFIELD brand skull pin is used with a MAYFIELD Skull Clamp.

## ⚠ CAUTION! ⚠

CAUTION: 60 lbs. (267 N) of force is recommended for adults. **40 lbs (178 N) of force is recommended for children.** All force settings are at the user's discretion. **Never exceed 80 lbs. (356 N) of force.**

### Skull Clamp Patient Fixation

## ⚠ WARNING! ⚠

Failure to properly position patient and to fully tighten and secure all adjustable portions of this or any similar device may result in skull pin slippage and serious patient injury, such as scalp lacerations, skull fracture, or even death.

1. Unlock the Index Locking Knob so the rocker arm can rotate freely 360 degrees. The Index Locking Knob is unlocked by rotating the knob to align the indicator with the open lock symbol; see Figure 4.


Figure 4 Unlocking Index Locking Knob

2. Rotate the rocker arm into the desired position for placement of the skull pins; normally perpendicular to the floor for prone or supine positions or parallel to the floor for patients in the sitting position.

Note: The Torque Screw force gauge provides the surgeon with a visual display of the clamping force value. (Each line on the gauge is 20 lbs or 89 N).

This will also align the rocker arm parallel with the direction of the surgical force (See Table 1 Skull Clamp Position Diagram).

3. Depress both of the Release Buttons on the body-portion of the Skull Clamp (See Figure 5 below) to disengage the Ratchet Extension and slide the two halves of the clamp apart to allow placement around the patient's head.


Figure 5 Disengaging Ratchet Extension

4. Position the 2-pin rocker arm so the skull pins are equidistant from a center line of the head.  
(see Figure 6).
5. Seat the skull pins securely into the cranial table. Skull pins should enter the scalp at a 90-degree angle, perpendicular to the surface of the head.


Figure 6 Positioning Rocker Arm Skull Pins

**⚠ WARNING! ⚠**

Make sure the trajectory of the single-pin in the Skull Clamp intersects the center of the patient's head and the rocker arm to correctly secure the Skull Clamp to the patient's head.

**⚠ WARNING! ⚠**

Avoid the areas of the frontal sinus, temporal fossa, major blood vessels, nerves, previously restored or abnormally thin bone.

6. Increase Skull Clamp force by advancing the Torque Screw to a desired setting, indicated on the visual readout, then back off one quarter turn.

Note: The Torque Screw provides the surgeon with a visual display of the clamping force applied. Each scribed line indicates a 20 lbs. (89 N) increment of force. Refer to Table 2 below:

| Torque Screw Line Indicator | Clamping Force  |
|-----------------------------|-----------------|
| 1                           | 20 lbs. (89 N)  |
| 2                           | 40 lbs. (178 N) |
| 3                           | 60 lbs. (267 N) |
| 4                           | 80 lbs. (356 N) |

Table 2 Torque Screw Line Indicator/Clamping Force


Figure 7 Standard Torque Screw

Note: To achieve independent fixation for adults, it is recommended that a minimum of 60 lbs. (267 N) of force be applied by the Skull Clamp. An additional force range reaching 80 lbs. (356 N) is available if desired.

7. Once the force is applied to the clamp and the patient is positioned for surgery, the Index Locking Knob should be turned to the locked position shown in Figure 4.

**⚠ CAUTION! ⚠**

Never exceed 80 lbs. (356 N) of force.

**⚠ CAUTION! ⚠**

Over-tightening the Index Locking Knob can result in damage to the skull clamp.

## **Attachment of The Skull Clamp to The Base Unit**

Once the Skull Clamp is applied to the patient's skull, the surgeon will maneuver the patient to the surgical position that is required for the procedure. With the patient in this position, the surgeon will hold the patient's head and the Skull Clamp and request that the components for the Base Unit to be brought up for attachment.

1. The MAYFIELD 2 Base Unit (A3100, A3101, A3102), should be attached to the operating room table with the Locking Handle unlocked and the Swivel Adaptor loose. In the event that other MAYFIELD Base Units (A2009, A2101, A2102) are going to be used with the Skull Clamp, be sure to also assemble a Swivel Adaptor (A1018, A2008) to complete the support system for the Skull Clamp. Make certain to securely attach the base unit to the operating room table.
2. The attached but loose components of the Base Unit should be brought up for attachment to the Skull Clamp.
  - Insert the mounting screw of the large starburst of the Swivel Adaptor into the large starburst of the Skull Clamp.
  - Turn the handle clockwise and tighten.
  - Care should be taken to maintain the position of the patient's head as requested by the surgeon.
3. After securely fastening the Skull Clamp to the Swivel Adaptor turn the Index Locking Knob to the locked position. Continue by checking all the locking knobs and locking handles of the Adaptor and Base Unit making certain all starburst teeth are meshed (where applicable) on all joints of the Base Unit after adjustments are complete.

## ⚠ CAUTION! ⚠

Before fully tightening, always be certain that the starburst teeth of the Swivel Adaptor and the Transitional Member fittings are the same size and properly mesh. Failure to do so may damage fittings. Figure 8 shows a typical starburst connection and proper meshing of teeth.


Figure 8 Mounting Skull Clamp

### **Maintenance, Care and Service Notes:**

#### **Maintenance and Care:**

To ensure proper function and to extend the life and performance of the equipment, Integra LifeSciences recommends the following:

| Recommended Action                                                                                                 | Recommended Frequency |
|--------------------------------------------------------------------------------------------------------------------|-----------------------|
| Return the device to the Integra LifeSciences Service and Repairs Department for detailed inspection and servicing | Once/ year            |
| Request that Integra Neuro-specialists perform routine inspections of the device                                   | Twice/ year           |

Note: In the absence of proper care and servicing of the device, negative effects may be seen after repeated reprocessing over time which may lead to reduced performance.

**Contact Information:** See the Service and Repair section for contact information on how to return your device.

**NOTE:** Any serious incident that has occurred in relation to the device for the user and/or the patient should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

**Device Disposal:** Follow hospital procedures for disposal of this device.

### **Service Notes:**

A routine equipment inspection made well before each case will help avoid problems. This check should include the following:

1. Store the Skull Clamp with the Index Locking Knob in the UNLOCKED position. This practice will provide the best performance from the knob.
2. While rotating the rocker arm through 360 degrees, lock and unlock the Index Locking Knob. To ensure proper locking engagement, the Index Locking knob must always be turned to the lock position. If the Index Locking Knob stops short of the locked position, a slight additional rotation of the rocker arm should complete the internal alignment and allow the locking mechanism to engage. If difficulties are encountered, the unit should be returned for service.
3. Check the Ratchet Extension release buttons to be sure that they allow for the Extension to disengage freely. Press them and move the Extension out of the body portion of the clamp - do this several times to make sure that it works smoothly and easily.

## Cleaning and Disinfection of the MAYFIELD 2 Skull Clamp

### Cleaning:

#### Option 1: Manual Cleaning

#### **CAUTION!**

Care should be taken to read completely this section on Cleaning, Disinfection and Decontamination of your MAYFIELD 2 equipment. It is important to know what method of cleaning/disinfection is needed based on what type of debris exposure the equipment has received. This is especially important when there is the possibility of exposure to a persistent pathogen that is resistant to normal cleaning and disinfection methods. Equipment with this exposure should be processed per the steps in the Decontamination Section. Integra does not make any claims regarding the effectiveness of the decontamination processes listed in deactivating pathogens, but rather, we are indicating that the device can withstand these procedures with minimal loss of function.

IMPORTANT: If the product has been exposed to persistent/ highly contagious pathogens<sup>1</sup>, DO NOT apply glutaraldehyde, alcohol or formalin to the device until the decontamination process is complete. (See Decontamination Section below for instructions)

The following steps are recommended for Manual Reprocessing of the MAYFIELD 2 Skull Clamp (A3059).

1. Disassembly: Separate the Ratchet Extension from the body portion of the Skull Clamp.

#### **WARNING!**

Visually confirm that there are no Skull Pins in the rocker arm or Torque Screw. It is recommended to completely remove the ratchet extension from the clamp before removing Skull Pins.

- Remove the Torque Screw from the Ratchet Extension.
2. Soak each component in a warm enzymatic solution (35°C–45°C) prepared in accordance with the manufacturer's instructions for 30 minutes.
  3. Wash and Rinse:
    - Scrub each component with a soft brush. The Skull Clamp can withstand detergents that fall within pH ranges of 3–11.
    - Clean thoroughly to remove any traces of blood and/ or debris. This is to prevent such blood and/ or debris from interfering with function or movement.
    - Inspect all components to ensure each component is visibly clean. Repeat as needed.
    - Rinse thoroughly with clean water (e.g. deionized water) to remove all traces of detergent.

<sup>1</sup> Persistent pathogens are resistant to deactivation by normal disinfection procedures and may require more aggressive treatment to achieve disinfection. See guidelines in Decontamination section for recommended procedures.

**⚠ WARNING! ⚠**

Use of highly alkaline detergents/ agents will cause a material degradation if present during autoclaving. Ensure device is thoroughly rinsed and all detergents are removed prior to use.

4. Dry all components with a clean, soft towel.
5. Re-assembly: After components are completely dry, re-assemble the components, leaving the Skull Clamp ready for the next use.
  - The Index Knob should be in the UNLOCKED position.
  - The Torque Screw should be inserted with the first thread showing.
  - The Ratchet Extension should be inserted into the body portion of the Skull Clamp.

**Option 2: Automatic Wash Cycle****⚠ CAUTION! ⚠**

Any deviation from this guideline could result in damage to equipment as well as improper cleaning results.

1. Disassemble the device as indicated in the Option 1 Manual Cleaning section of the MAYFIELD 2 Skull Clamp.
2. Rinse components under warm tap water prior to placing in washer.
3. Arrange device in a way to avoid contact between components.
4. Select instrument cycle and ensure proper programming. (Reference Table 3)

| Phase                  | Time in minutes | Water Temperature         | Detergent type & Concentration                                                        |
|------------------------|-----------------|---------------------------|---------------------------------------------------------------------------------------|
| Pre-Wash 1             | 04:00           | Cold tap H <sub>2</sub> O | N/A                                                                                   |
| Enzyme Wash            | 04:00           | Hot tap H <sub>2</sub> O  | Neutral pH enzymatic detergent (e.g. Endozyme® AW Triple plus w/ A.P.A.), 1:128 ratio |
| Wash 1                 | 10:00           | 60° C                     | Neutral pH cleaner, concentrated (e.g. Renu-Klenz™) 1:256 ratio                       |
| Rinse 1                | 00:30           | Hot tap H <sub>2</sub> O* | N/A                                                                                   |
| <i>Thermal Rinse**</i> | 02:00           | 82.2° C                   | N/A                                                                                   |

Table 3 Instrument Cycle

\* No limit on tap water temperature.

\*\* Optional phase for disinfection of components - minimum water temperature as indicated or per washer manufacturer specifications for the thermal rinse cycle.

5. Remove from washer and dry with a clean, soft towel.
6. Inspect all components to ensure there is no visible organic debris or residue from cleaning agent. Repeat process if any visible soil is detected.

### Additionally Validated Cleaning Cycles

| (Pre-Wash)                                                      | Wash        |                |                                                                   | (Rinsing)                             | Thermal Disinfection | (Drying) (Optional) |
|-----------------------------------------------------------------|-------------|----------------|-------------------------------------------------------------------|---------------------------------------|----------------------|---------------------|
|                                                                 | Temperature | Hold Time      | Detergent Treatment                                               |                                       |                      |                     |
| Cold water from building supply (less than 110° Or 43°C)/4 min. | 93°C        | 1 or more min. | Alkaline (e.g. neodisher® SeptoClean 1%)                          | Hot water from building supply/2 min. | Optional             | 90°C/20 min.        |
| 50°C/7 min                                                      | 90°C        | 1 or more min. | Neutral Enzymatic (e.g. Enzol® per manufacturer's recommendation) | Hot water from building supply/2 min. | Optional             | 90°C/20 min.        |

Table 4 Additionally Validated Cleaning Cycles

### Optional Thermal Rinse:

- A Thermal Disinfection rinse may be added after the first rinse cycle as indicated in Table 3.

**Disinfection:**

The MAYFIELD 2 Skull Clamp must be thoroughly cleaned after each use. Where possible, keep the soiled areas moist until the disinfection process can be initiated. Disinfection of device components following normal usage may be achieved using the following methods:

**Option 1: Thermal Rinse:**

A Thermal Disinfection rinse may be added after the first rinse cycle as indicated in Table 3. See Cleaning section.

**Option 2: Chemical / Bleach Soak:**

 **WARNING!** 

Bleach solutions produce chlorine gas. Prepare in a well-ventilated area.

Note: Hypochlorite potency decreases as bleach ages. Use only recently purchased/ opened bleach to prepare disinfection solutions. Do not use pre-diluted bleach. Prepare a new bleach solution for each application.

Hypochlorite (Bleach) mixing guidelines to obtain 20,000 ppm solution:

| Bleach concentration | Preparation Instructions                     |
|----------------------|----------------------------------------------|
| 4.5%                 | Pour 4 units of bleach into 5 units of water |
| 5%                   | Pour 2 units of bleach into 3 units of water |
| 6%                   | Pour 1 unit of bleach into 2 units of water  |

Table 5 Hypochlorite Bleach Mixing Guidelines

1. Prepare a solution for general disinfection (See Table 5 and Table 6 guidelines for bleach solution)

| For General Disinfection                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"><li>• Verify Concentration of base hypochlorite is 4.5% or more</li><li>• Mix ½ unit of bleach base solution into 10 units of water</li></ul> |

Table 6 Hypochlorite Bleach Mixing Guidelines

2. Disassemble the device. (see above disassembly instructions in Option 1 Manual Cleaning)
3. With the device disassembled, place components in a shallow pan and add enough bleach solution to cover the device completely.
4. Soak components for 15 minutes.
5. Rinse thoroughly with clean water. Ensure no bleach residue remains on components.

 **CAUTION!** 

Use of highly alkaline detergents/ agents will cause a material degradation if present during autoclaving. Ensure device is thoroughly rinsed and all detergents are removed prior to use.

6. Dry with a soft, clean towel.

### Option 3: Steam Sterilization/Autoclave:

#### Prepare for Sterilization:

1. Visual Inspection:
  - Visually inspect the Skull Clamp components for cleanliness.
  - Ensure all components are free of organic debris and are in proper working order. If damage is observed on any components, the MAYFIELD 2 component should be returned immediately to the Integra Service and Repair Center for servicing. (Please refer to the Maintenance and Servicing Section.)
2. Wrap Components for Sterilization:
  - Ensure the components are completely dry prior to wrapping for sterilization. Components may be dried with a clean, soft towel or cloth.
  - Individually wrap each component to be sterilized, (partially disassemble per the above disassembly instructions for cleaning) and swivel adapter or Tri-Star adapter using an approved sterilization wrap appropriate for a steam sterilization/ moist heat sterilization method.

#### Sterilization:

1. Use one of the below described sterilization cycle parameters. Place pre-wrapped components in a validated steam sterilizer. (Reference Table 7)
2. Follow the sterilizer manufacturer's instructions for loading and unloading the devices and operating the sterilizer.

Note: Equipment cycles may vary due to sterilizer size and maximum load capacity. To ensure sterilization, use the sterilizer manufacturer recommendations in regards to maximum capacity to validate individual equipment cycles.

#### **CAUTION!**

Avoid overcrowding on sterilizer shelves/ racks and heavy loads. Reference the sterilizer manufacturer instructions for the proper use while operating the sterilizer.

#### **CAUTION!**

Prolonged exposure time may affect product life. Carefully inspect for any visible damage or malfunction after sterilization.

#### **WARNING!**

Ensure device has properly cooled prior to removing from the sterilizer. Follow the sterilizer and established procedures when operating the equipment.

The MAYFIELD 2 Skull Clamp (A3059) has been validated for the following sterilization cycle parameters:

### Steam Sterilization Validated Cycles

| Cycle Conditions                            | Sterilization Temperature        | Exposure Time (minutes) | Dry Time (minutes) |
|---------------------------------------------|----------------------------------|-------------------------|--------------------|
| Wrapped, Pre-vacuum<br>US cycle             | 132°C -0/+3°C<br>(270°F -0/+5°F) | 3                       | 30                 |
| Wrapped, Pre-vacuum<br>US- Alternate Cycle  | 132°C -0/+3°C<br>(270°F -0/+5°F) | 4                       | 60                 |
| Wrapped, Pre-vacuum<br>UK/EU Cycle          | 134°C -0/+3°C<br>(274°F -0/+5°F) | 3                       | 30                 |
| Wrapped, Pre-vacuum<br>German Cycle         | 134°C -0/+3°C<br>(274°F -0/+5°F) | 5*                      | 60                 |
| <i>Wrapped, Pre-vacuum<br/>French Cycle</i> | 134°C -0/+3°C<br>(274°F -0/+5°F) | 18*                     | 60                 |

Table 7 Validated Cycles

Note: It is the responsibility of the end user to validate alternate sterilization methods and cycles that do not comply with the validated sterilization methods and cycles specified within this manual.

### **Decontamination after Exposure to Highly Contagious Diseases (WHO):**

When the MAYFIELD 2 Skull Clamp has been exposed to pathogens that resist normal disinfection methods, the product should be handled more carefully:

- Wear appropriate personal protective equipment (eye/face protection, gloves)
- Keep soiled areas of the device moist. Do not apply gluteraldehyde, alcohol or formalin to the device until the decontamination process is complete.

Note: Do not subject the Skull Clamp to automated cleaning processes until after the decontamination procedure.

1. Remove Ratchet Extension from the Clamp Base (see Figure 5).
  - Visually confirm there are no skull pins in the Rocker Arm or Torque Screw sockets.
2. Remove and dispose of skull pins if present.
3. Remove gross debris and soil from the clamp components using soft cloth or paper towels. Minimize formation of aerosols or droplets during the initial preparation of the product for decontamination.
4. Remove Torque Screw from Ratchet Extension.

---

\* Steam Sterilization is validated between 132°C - 137°C with an exposure time between 4 - 18 minutes and a dry time of 60 minutes.

5. Place parts in a shallow pan and add enough chlorine solution to cover the device completely.

- Allow to soak for 60 minutes.

6. Hypochlorite/ Bleach Soak:

### **CAUTION!**

Use of highly alkaline detergents/ agents will cause a material degradation if present during autoclaving. Ensure device is thoroughly rinsed and all detergents are removed prior to use.

- Prepare a hypochlorite/bleach solution. (Refer to Chemical/ Bleach soak instructions for mixing instructions)
- Place parts in a shallow pan and add enough bleach solution to cover device completely.
- Soak for 60 minutes
- Remove the device from the bleach solution and rinse components completely with water.
- Dry all components with a clean, soft towel

7. Sterilization/Autoclave

- The Skull Clamp may be autoclaved at 134° C for 18 minutes to 1 hour. The Skull Clamp has been validated to withstand this process up to 15 times without loss of performance.

Note: Integra makes no claims that the following protocol is effective in neutralizing any specific pathogens.

#### **World Health Organization (WHO) Protocol for Decontamination of highly infectious diseases<sup>2</sup>:**

“Immerse [instrument] in sodium hypochlorite for 1 hour; remove and rinse with water, then transfer to open pan and heat in a gravity displacement autoclave at (121° C) or porous load (134°C) for 1 hour.”

The life expectancy of the MAYFIELD products is expected to be as long as 7 years.

<sup>2</sup> WHO/CDS/CSR/APH/2000.3, “WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies”, (Switzerland, March 1999) pp 29

## **Integra Standard Warranty**

INTEGRA warrants to the original purchaser only that each new MAYFIELD product is free from manufacturing defects in material and workmanship under normal use and service for a period of one year (except as otherwise expressly provided as to accessory items) from the date of delivery by INTEGRA to the first purchaser, but in no event beyond the expiration date stated on any product labeling.

- Surgical instruments are guaranteed to be free from defects in material and workmanship when maintained and cleaned properly and used normally for their intended purpose.
- Any covered product that is placed by INTEGRA under a lease, rental or installment purchase agreement and that requires repair service during the term of such placement agreement shall be repaired in accordance with the terms of such agreement.

If any covered defect occurs during the warranty period or term of such placement agreement, the purchaser should communicate directly with INTEGRA's home office. If purchaser seeks to invoke the terms of this warranty, the product must be returned to INTEGRA at its home office. The defective product should be returned promptly, properly packaged and postage prepaid. Loss or damage in return shipment to INTEGRA shall be at CUSTOMER's risk. INTEGRA's sole responsibility under this warranty shall be repair or replacement, at INTEGRA's sole discretion at INTEGRA's expense, subject to the terms of this warranty and applicable agreements.

IN NO EVENT SHALL INTEGRA BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA shall not be responsible for, any loss arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA. THIS LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS.

INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS. No warranty or guarantee may be created by any act or statement nor may this Standard Warranty be modified in any way, except as a result of a writing signed by an officer of INTEGRA. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct.

## **Service and Repair**

For service and repairs within the United States, send all instruments for service or repair to:

Integra LifeSciences Corporation  
4900 Charlemar Drive, Building A  
Cincinnati, OH 45227, USA

(Always include the purchase order number and a written description of the problem)  
Or phone: 877-444-1114

For service or repair outside the United States, contact your local authorized Integra representative or customer service:

International: +33 (0) 4 37 47 59 50 · +33 (0) 4 37 47 59 25 (fax) · csemea@integralife.com

France: +33 (0) 4 37 47 59 10 · +33 (0) 4 37 47 59 29 (fax) · custsvcfrance@integralife.com

Benelux: +32 (0) 2 257 41 30 · +32 (0) 2 253 24 66 (fax) · custsvcbenelux@integralife.com

Germany: +49 (0) 2102 5535 6200 · +49 (0) 2102 553 6636 (fax) · custsvcgermany@integralife.com

United Kingdom: +44 (0) 1264 312 725 · +44 (0) 1264 312 821 (fax) · custsvcsuk@integralife.com

Switzerland: +41 (0) 2 27 21 23 00 · +41 (0) 2 27 21 23 99 (fax) · custsvcsuisse@integralife.com

Italy: +39 02 577 8921 · +39 02 575 11 371 (fax) · custsvcitaly@integralife.com

U.S. Patent 9,681,924

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation in the United States and/or other countries. MAYFIELD is a registered trademark of SM USA, Inc. and is used by Integra under license. Endozyme is a registered trademark of Ruhof Corporation. Renu-Klenz is a trademark of Steris Corporation. neodisher is a registered trademark of Chemische Fabrik Dr. Weigert GmbH & Co. Enzol is a registered trademark of Johnson & Johnson.

©2021 Integra LifeSciences Corporation. All rights reserved. 451A3059 Rev FAB 2021-06 0651223-3

This page intentionally left blank.

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Integra LifeSciences Corporation  
4900 Charlemar Drive, Building A  
Cincinnati, OH 45227, VS

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Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation in the United States and/or other countries. MAYFIELD is a registered trademark of SM USA, Inc. and is used by Integra under license. Endozime is a registered trademark of Ruhof Corporation. Renu-Klenz is a trademark of Steris Corporation. neodisher is a registered trademark of Chemische Fabrik Dr. Weigert GmbH & Co. Enzol is a registered trademark of Johnson & Johnson.  
©2021 Integra LifeSciences Corporation. All rights reserved. 451A3059 Rev FAB 2021-06 0651223-3

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©2021 Integra LifeSciences Corporation. Alle rechten voorbehouden. Publicatienummer 451A3059 Rev FAB 2021-06 0651223-3